DR
**SEBBAH
SIMON**

Aesthetics doctor in Paris
Consultant for
Doctor and Co medical centre
Former house doctor
for the Paris Hospitals

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ATTIVA ENDODERMAL RADIOFREQUENCY

AN INNOVATION FOR FIRING THE SKIN ON THE
FACE AND BODY S

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External RF is well known and commonly used in treatments designed to firm up the skin, reduce wrinkles and reduce cellulite.

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This newest innovation from Attiva uses a technology that delivers internal radiofrequency subcutaneously, allowing us to go even further in the skin firming process because we can target the deep skin tissue more precisely.

It allows us to treat sagging skin on the face (oval, submental region) but also sagging skin on other parts of the body (arms, around the bellybutton after significant weight loss or after pregnancy, thighs, above the knees). This internal radiofrequency technology is emitted through the tip of a cannula that is inserted underneath the skin. The Attiva radiofrequency system heats up the sagging tissues to a temperature of between 48 and 52°C, allowing for the stimulation and synthesis of collagen and elastin, which improves the skin's firmness. //

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SAFETY AND EFFICACY:

This technology has a double temperature control system that enables it to be both efficient and, above all, not to cause any damage to the skin treated.

- The internal temperature is monitored and measured at the tip of the cannula, which allows us to stay within the most efficient temperature window (48° to 52°).
- The external temperature is monitored and measured in real time using a thermal camera placed above the treatment area.

This allows us to ensure that the external temperature of the skin always stays below 42°, thus avoiding any risk of burning the skin.

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“ITS VERSATILITY IS VERY INTERESTING, AS IT CAN BE USED PURE AS A VOLUMISER OR DILUTED TO FIRM AND DENSIFY THE SKIN OF THE ARMS, THIGHS AND BUTTOCKS.”

PROTOCOLE :

The treatment is pain-free because it is carried out under a light local anaesthetic, and it only takes from 45 minutes to an hour per treatment zone.

The advantage of this treatment is that there are little or no side effects for the patient, and they can get back to work very quickly.

When the patient leaves the appointment, there may be some slight swelling in the treatment area, which is mainly due to the anaesthetic injected prior to treatment.

The patient might feel some sensitivity in the treatment area for 2 to 3 days afterwards.

RESULTS:

The results gradually become visible in the weeks that follow. A check-up is scheduled 3-4 months after treatment and another session may be recommended at this time, depending on the initial condition of the sagging and our desire to firm up the treatment zone further. This makes the results last longer. A top-up session might be scheduled every 2-3 years. ●